Sebond Dentin Bond Opaker

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Sebond Dentin Bond Opaker

Intended use:

Framework covering materials are light-curing or chemically curing composite-based materials for covering different framework materials in preparation for subsequent veneering with veneering composite or denture resin.

Product description:

Sebond Dentin Bond Opaker is a light-curing, tooth-colored paste opaque with an integrated adhesive component for coating frame- work retentions and for clasp coating.

General product description:

– Light-curing one-component material for color-coating frame- work retentions and clasps made from commercially available dental metals and alloys. The material is light-curing and is to be polymerized in a light polymerization unit. Sebond Dentin Bond Opaker contains an additional adhesive component which assures the bond to the framework material. The use of an additional bonding system is therefore not necessary.
– The material is applied to retentions on frameworks before their completion with acrylic to coat them in tooth color.

Composition:

– Silicon dioxide, glass powder, di-urethane dimethacrylate, 2-[[[(butylamino)carbonyl]oxy]ethyl acrylate, tetramethylene dimethacrylate, pigments, initiators
– Filler content: 45 % by weight; anorganic filler (0.005 - 0.3 µm)

Indications:

For coating framework retentions (commercially availabe dental metals/ alloys) and clasp coating for new restorations and repairs.

Contraindications:

– Do not use in case of a known allergy against on of the compo- nents.
– Sebond Dentin Bond Opaker is not suited to composite veneering.

Interactions with other products:

Phenolic substances (such as e. g. Eugenol) inhibit polymerization. Do not use materials which contain such substances.

Application / processing:

Processing time of the materials: 1-3 minutes, depending on the light conditions.

Preparatory tasks::

– **After casting**, the clasp or framework is designed, cast, milled with CAD/CAM and finished as usual, according to the applicable guidelines for dental technology.
– Suitable framework materials are commercially available dental metals and alloys, PMMA, PEEK and zirconium dioxide. All materials shall be processed and prepared for further processing in accordance with the manufacturer's instructions.
– The surface is sand-blasted with 50 µm aluminium oxide (2 bar) and is blown with oil-free compressed air.

Procedure:

Sebond Dentin Bond Opaker is applied thinly (wash-like) with a firm short-hair brush. At least two layers must be applied to obtain a complete coating. An insufficiently covering opaque layer negatively influences the color effect. The flowability of the opaque can be optimized by stirring on a mixing block. Make sure to apply the material in an especially thin coat around the retentions. Poly- merize each layer separately (see table "polymerization times"). The opaques can be intermixed for an individual color effect.

Subsequent tasks:

Finally, the acrylic for covering the retentions can be applied. In case of a clasp coating, finish the restoration as usual.

Troubleshooting / FAQ list

– Clean the surfaces of the framework from grease and polishing residue.
– Do not mix Sebond Dentin Bond Opaker with opaque liquids or opaque powders.
– The opaque must only be applied in very thin layers.
– Opaque layers which are too thick prevent an optimal polymer- ization and thus weaken the bond. Extending the polymerization time will result in a better curing depth.
– Do not use intermediate polymerization systems for the poly- merization of this product.

Error	Cause	Remedy
does not cure	layer too thick	only apply wash- like in several layers
	insufficient polymerization	observe poly- merization times – check lamp, replace if necessary – extending the polymerization time will result in a better curing depth – do not use an intermediate polymerization system
opaque chips	opaque not polymerized sufficiently/ applied too thickly	nonly apply wash- like respectively in very thin layers
		incorrect poly- merization times/ check lamp/ replace lamp if necessary
material in the syringe is too solid	syringe stored at a temperature below 10 °C	store material at room temperature or let the material reach room tem- perature in time for its use
greasy surface	insufficient polymerization	observe poly- merization times
		check system/ do regular maintenance

General information:

Use a suitable light polymerization system with an emission spectrum of at least 310 - 500 nm. The required physical properties are only achieved with perfect lamps. Regular checks of the light intensity according to the manufacturer's specifications are required. Do not use intermediate polymerization systems.

Storage:

– Store at 10 - 25° C (50 - 77° F).

– Close syringe well after each use.

Shelf-life:

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

Hazard statements / precautionary statements:

Contains tetramethylene dimethacrylate, diphenyl(2,4,6-trimethyl- benzoyl)phosphine oxide, di-urethane dimethacrylate, 2-[[[(buty- lamino)carbonyl]oxy]ethyl] acrylate. May cause allergic skin reactions. Avoid breathing dust/ smoke/ gas/ vapor/ steam/ aerosol. Wear protective gloves.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, the composition of the medical device is available upon request. The dentist should con- sider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the attending physician with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Polymerization times:

Light polymerization system	Time
Spekra LED	1 min.
Spekra 2000	3 min.
HiLite / Uni-XS	90 sec.
Labolight LV-II / III	1 min.
Solidilite	1 min.

Disposal:

Disposal of the content/ of the container in accordance with local/ regional/ national/ international regulations.

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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